

# ENGINEERED FOR PRECISION DESIGNED FOR CONFIDENCE



HIGH PERFORMANCE SURGICAL CLIPPER WITH NON-CONTACT  
BLADE REMOVAL FOR ENHANCED HYGIENE AND ZERO  
CONTAMINATION RISK

Pre-op hair removal with speed, safety,  
and hygiene.



**180 mins runtime, 120 mins charge**



**IPX7 water-resistant**



**Non-contact blade ejection**



**Universal blade for all applications**



**Precision blade design minimizes  
skin contact reducing trauma**



**Digital LED battery status level indicator**



**UK Design &  
Development**

## Why Choose Our Surgical Clipper?



**Unrivaled Motor Power** : High torque for fast, efficient trimming



**One Universal Blade** : Precision cutting across body areas



**Ergonomic Grip** : Better handling, less risk of slips



**Whisper Quiet (65dB)** : Comfortable experience for patients



**Water Resistant** : Easy to clean and disinfect



**Non-Contact Blade Ejection** : Reduces cross-contamination




**Powerful Lithium Battery** : 3 hrs runtime, quick charge

## Professional Recommendation

Worldwide all infection control guidelines, including those from the CDC, NICE, WHO, AORN, and Australian health authorities, recommend electric clippers with single-use blades as the optimal technique for preoperative hair removal.

- Manual shaving causes skin damage, increasing infection risk.
- Clipping is preferred over shaving as it causes less skin trauma and is associated with lower surgical site infection (SSI) rates.
- Timing hair removal immediately before surgery further reduces infection risks.
- Proper training of staff in clipping technique and hygiene is critical for optimal outcomes.



## Main Benefits for CLINICIANS & H.C.W

### Reduces Risk of Cuts and Micro-Injuries

Safe blade design prevents direct skin contact.

### Quick, Simple Operation :

Streamlines patient preparation and saves crucial time.

### Clear Incision Sites :

Facilitates better dressing adhesion and reduces infection rates.

### Reliable Performance:

High-efficiency motor and battery ensure readiness procedure after procedure.

---

## Main Benefits for Patients

### Minimal Skin Irritation :

Low risk of micro-injuries, infection, and irritation.

### Comfortable Experience :

Low noise and an efficient, pain-free process.


### Safe for Delicate Areas :

Specially designed blade and handle for use in sensitive regions.

---

## Simple Steps for Use

- 1 Insert blade until it clicks
- 2 Hold like a pen and glide against hair
- 3 Use eject buttons for blade removal
- 4 Rinse, disinfect, and charge for next use



# What's Included in the Box

Surgical  
Clipper Unit

Charging  
Adapter  
(100-240V)



Instruction  
Manual

## ORDERING INFORMATION

REF Order number	Product name	Units per Master Carton	Master Carton
LP003	LP Surgical Clipper Handle with Charging stand & Adapter	20 nos	10 kgs (0.040 cbm)
LP004	LP Surgical Clipper Blade	6 x 50 unit ( 300 nos )	5.2 kgs (0.049 cbm)

### Contact & Address :

**L. P. Corporation Ltd**

71-75 Shelton Street, Covent Garden, London WC2H 9JQ,  
UNITED KINGDOM

✉ info@lpcorporation.com

🌐 www.lpcorporation.com



**Delivering Innovation. Shaping Hygiene**